SELEÇÃO

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- B0111-100 (P) - 408 x 408 x 630 mm
- B0111-210 (P) - 408 x 408 x 1285 mm
- B0111-300 (P) - 605 x 605 x 880 mm
- B0111-500 (P) - 605 x 605 x 1285 mm

Portas

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Aplicação

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Sistema de Vácuo

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
- Commonly used in CSSD

Vacuum System

- Vacuum pump – one stage – 1,5 hp (available for all volumes)
- Vacuum pump – one stage – 3,0 hp (unavailable for 100 l)

Steam Source

- External Source (See standard requirements)
- Electric steam generator
- Switch (direct steam + steam generator)

Panels

- Removable side panel
- Back panel (single door)
- Top panel

Optional components

- Barrier panel
- Residual current circuit breaker
- Built-in air compressor

Electrical power

- 220V 380V
- 60Hz 50Hz

Accessories

- Load rack
- Rack trays (available for 100 l volume only)
- Wall support
- Preparation bench
- Rack trolleys
- ISO baskets
- Floor support

Installation Components

- Water treatment Air compressor
- Air treatment Sanitary filter
- Steam treatment Softener
- Osmosis

Services

- Installation Start Up
- IQ (Instal. Qual.) OQ (Oper Qual.)
- PQ (Perf. Qual.) Training
- Maintenance Parts kit

FOTO ILLUSTRATIVA - DETALHES MODEM

APPLICATION

Sterilizers HI SPEED II line using steam under pressure at high temperatures as a sterilizing agent in cycles by pulsing air removal and automated vacuum control system through dedicated type microprocessor controller PLC industrial. Built with the use of the best technology available in the areas of mechanics, thermodynamics, hydraulics, electronics and software, equipment is robust, secure and high performance. On models with one or two doors and with capacities ranging from 100 to 500 liters are intended for material sterilization sector in hospitals and other healthcare facilities, medical laboratories, quality control laboratories and production processes of pharmaceutical and others. Its use is suitable for sterilization of packaged porous material, instruments and utensils packaged or not, glassware, gloves, syringes and other materials in rubber or alloys resistant to high temperatures.

MAIN FEATURES

High performance: Equipment designed and tested according to the NBR ISO 17665 standard, ensures reliable, accurate and valid sterilization processes.

Speed and Security: Ensured by precise control routines implemented by proprietary software and high precision instrumentation.

Simplified operation: Microprocessor control and intuitive interface screen with "Touch Screen" colorful. Facilitates adherence of users to process requirements.

Versatility: Flexibility in scheduling processes. Access to programming and equipment functions through user passwords with configurable permission levels.

Sustainable Processes: Environmentally clean

and safe equipment without residues to cycle completion. Chambers, doors and vapor transport pipes with thermal insulation. Energy saving, lower heat dissipation and lower thermal load on the air conditioning system.

Installation and Maintenance: Facilitated by panel for electrical components, rear and one side.

MODELS AND CAPABILITIES

The HI SPEED II family consists of small, medium and large sterilizers, with the capacities below:

- B0111-100 - Capacity of 100 l
- B0111-210 - Capacity of 210 l
- B0111-300 - 300 l capacity
- B0111-500 - Capacity of 500 l

STANDARDS

- **EN ISO 17665-1** - Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.
- **NBR 11816** - Vacuum steam sterilizers for healthcare products.
- **ASME, Section VIII, Division I** - Boiler and pressure vessel code.
- **NR 13** - Boilers and Pressure Vessels.

CONSTRUCTION

Sterilization Chamber - Built in double wall (inner chamber and outer chamber). Inner chamber with circular profile built in AISI 316L stainless steel, with a thickness of 3/16" (4,76 mm) and sanitary polished finish. Outer chamber in circular profile with the same material as the inner chamber. The set is sized according to ASME standard to withstand the following pressures:

- Maximum working pressure of 3,0 kgf/cm².

- Hydrostatic test pressure of 4,5 kgf/cm² - 1,5 times the PMTA.

A 1" diameter validation input is available on the side of the equipment enabling the introduction of temperature and pressure sensors in operational and performance qualification processes.

Externally, a chloride-free mineral fiber layer reduces vapor condensation and heat radiation. Inner chamber drain protected with 1" diameter stainless steel perforated plate filter located near the load side face.

Alternatives for the Chamber Material

- AISI 316 L - 3/16" (4,76 mm)

Finishes

- Polished Sanitary

Structure: The chamber assembly is mounted on a carbon steel structure with anti-corrosive paint and adjustable feet height to allow for perfect leveling of the equipment.

Pipes: the pipes welded in the chamber are constructed of AISI 316 stainless steel.

The pipes on the chamber casing are made of brass. Connections are threaded.

Alternatives to Pipes

Outside coating

- Brass/ Bronze

The chamber assembly is protected by a safety valve built in bronze, sealed and with a verification mechanism, with a threaded connection.

Alternatives for Safety Valve

- Cast Iron/ Bronze

Panels and Barrier panel

Front panel and removable side panel manufactured in AISI 304 stainless steel and brushed finish or carbon steel with anti-corrosive paint. As standard the equipment is sold without the top and rear finishes.

A panel (optional) fixed to the equipment chamber completes the sealing between the environments where the sterilizer will be anchored. This seal can meet hermetic requirements between controlled areas or just guarantee the maintenance of the pressure differential between contiguous areas, as required by hospitals.

Alternatives panels

Front panel

- Stainless steel AISI 304 brushed

- Stainless steel without nickel AISI 444 brushed

- Carbon steel with anticorrosion protection and painting

Side panel

- Removable stainless steel AISI 304 brushed

- Removable stainless steel without nickel AISI 444 brushed

- Removable carbon steel with anticorrosion protection and painting

Superior finish

- Brushed AISI 304 stainless steel

- Brushed Nickel-free AISI 444 stainless steel

- Carbon steel with anticorrosive protection and paint

Top panel

- Stainless steel AISI 304 brushed.

- Stainless steel without nickel AISI 444 brushed.

- Carbon steel with anticorrosion protection and painting.

Back panel (single door)

- Stainless steel AISI 304 brushed.

- Stainless steel without nickel AISI 444 brushed.

- Carbon steel with anticorrosion protection and painting.

Barrier - hospital (double door)

- Stainless steel AISI 304, screwed.

Doors: The equipment can be supplied with one or two doors when installed in a sanitary barrier, with automatic opening and closing, activated through the command interface. The door is vertically moved via a pair of pneumatic pistons, rising on closing and lowering behind the cabinet on opening.

Mechanical safety locks prevent the doors from opening during the sterilization process and the simultaneous opening of the loading and unloading doors in the case of equipment with two doors in the loading and unloading phases of the equipment. Sealing the door by a silicone gasket housed in a channel in the chamber, which at pressurized closing is closed against the internal face of the door and in the retracted opening by means of a vacuum pulse. Both the face of the door and that of the channel are rectify for proper adjustment.

The door(s) has an internal part constructed of AISI 316L stainless steel plate. Internal insulation

made with chloride-free mineral fiber mat, internal structure through steel profiles and external finishing by carbon steel plate with anti-corrosive protection or optionally in AISI 304 stainless steel with brushed finish.

Alternative for Door

- Single door

- Double door

Material (inner face)

- AISI 316L stainless steel

Finishing

- Polished Sanitary

Entry of clean air: to break the vacuum through hydrophobic sterilizable filter with 99,9997% efficiency for larger particles or equal to 0,22µm.

Vacuum System: For carrying out the vacuum pulses in the load conditioning at the beginning of the cycle and drying in the final stage of the process, through a single stage vacuum pump, monobloc type with water ring, sized to meet the requirements of current standards.

Options for Vacuum Systems

- Vacuum Pump - one stage

Vacuum pump power: depends on the volume of the chamber and the configuration options available for each model guarantee the successful completion of Bowie & Dick and hollow load tests as defined in national and international standards.

Available Vacuum pumps

Vacuum Pump - one stage

- 1,5 Vacuum Pump hp for 100, 210, 300 and 500 liter sterilizers.

- 3,0 Vacuum Pump hp for 210, 300 and 500 liter sterilizers.

The steam used in the sterilization of materials is preferably supplied through an electric steam generator incorporated into the equipment and fed by water with adequate quality for the process.

Steam supply - Alternatives

- With preparation for receiving electric steam generator.

- Without steam generator. Steam supplied from an external source.

NOTE: The generator must be set up as a separate item. See construction details of the generators in

“optional components”.

ELECTRIC PANEL

Control panel arranged in height that facilitates viewing and operation of the sterilizer load side.

Microprocessor controller dedicated PLC industrial type with routine self-tests, self-diagnosis and supervision cycles by “watchdog timer” (watchdog time) according to the current standard.

Human Machine Interface (HMI) Touch Screen with color screen 7”, IP65, allows you to monitor and control all machine functions and parameters of the cycle. With USB port for software update.

The Interface allows you to select and start a sterilization cycle, follow the steps of the process with the visualization of the temperature, pressure measurements, estimated time to end the phase; manage users and passwords with the possibility of registering 40 users with configurable access levels; configure end of cycle alert; view fault and security alarms; change system date and time; screen adjustment; access to the record of the last alarms with the possibility of exporting this data; and choice of operating language.

The interface also provides access to the operation routines for programming the cycle parameters in security bands defined in the project, access to manual operation, calibration and maintenance using a user login and password, and the programming of routine maintenance activities, depending on the number of cycles or hours of operation.

Secondary panel on the sterile side (discharge), in equipment with double door, equipped with a manometer with glycerin for reading the pressure of the internal chamber, operation button for the discharge door, emergency button and lamps for indicating the cycle in progress, loading door closed and active alarm.

Process control performed by the main PLC with temperature data sent by the PT-100 temperature sensor allocated next to the drain. Pressure data sent by the inner chamber and outer chamber transducers. The PLC has internal control of process times.

The entire sterilization process can be monitored in real time through the interface that informs the operator:

- Chamber temperature;
- Chamber pressure;
- Jacket pressure;
- Time elapsed from the exposure / sterilization phase;
- Time elapsed to complete the drying phase;
- Cycle status information;
- Real-time graphical display of temperature and pressure data.

The panel also contains the security and control instruments indicated below:

- On-off key;
- Emergency button;
- Manovacuumeter of the chamber;
- Jacket pressure gauge;
- Manometer of the steam generator or external steam;
- Printer for recording the parameters reached in the process

Instrumentation and control

Pressure control of the internal and external chambers performed by electronic pressure transducers made of stainless steel, with electronic temperature compensation and absolute pressure reading. Pressure indications are made digitally via the panel interface in bar.

Temperature control by PT-100 class A sensor, placed next to the chamber drain, it reads the process temperature. Temperature indications are performed digitally via the panel interface in °C.

Load sensor Sensor for measuring temperature in load, available as a control option, mainly in liquid sterilization cycles. With flexible cable for better positioning inside the load.

Paperless System for recording cycle information in electronic form (in PDF format) without the use of paper. The data can be accessed via USB and will only be recorded using a USB stick.

STERILIZATION CYCLES - PROFILE

The HI SPEED II Sterilizer Line can be supplied with three alternative and distinct programs with profiles with configurable cycles. In both, cycles with profile for density, liquids and lactations are available:

“Program A” is loaded as standard on sterilizers and provides the user with 8 sterilization cycles with density profiles from 101 to 135°C and 1 Bowie&Dick cycle.

“Program L” is loaded as standard on sterilizers and provides the user with 2 sterilization cycles with 121°C density profiles, 4 cycles with 110°C lacquer profiles and 3 cycles with 121°C liquid profiles.

“Program B” is loaded as standard on sterilizers and provides the user with 5 sterilization cycles with density profiles from 101°C to 135°C, 2 cycles with liquid profiles from 101°C to 135°C, 1 cycle from Bowie&Dick and 1 Leak-test cycle.

DENSITY MATERIAL CYCLES

Profile cycles for the sterilization of textile materials, glassware, metallic instruments packed in perforated boxes or containers, devices with lumen, thermo-sensitive materials in latex and / or plastics with resistance to sterilization temperature. The cycle consists of 06 distinct phases:

Phase 1 - Air removal and conditioning

Through alternating pulses of vacuum and steam. This phase is only configurable in the Bowie & Dick

cycle and this configuration is automatically loaded in other cycles with this profile. The parameters for this phase are:

- Initial vacuum pulse with a programming range between 0,05 and 0,90 bar absolute.
- Steam pulse with programming range between 1,0 and 2,0 bar absolute.
- Waiting time after steam pulse with a programming range between 01 and 9999 seconds.
- Number of pulses with a programming range between 1 and 9999.

Phase 2 - Load heating

Through steam injection until the sterilization temperature in the chamber and the load to be sterilized is reached. No programmable parameters. The equipment controls pressure and temperature in the chamber according to the sterilization temperature defined for the process.

Phase 3 - Exhibition

Includes a time interval of up to 60 seconds for the temperature to settle within the range defined for the process (sterilization temperature + 2 °C) and the threshold time with maintaining the process temperature in that same range until the end of the phase. The parameters for this phase are:

- Sterilization temperature with a programming range between 101 and 135 °C.
- Exposure time with a programming range between 1 and 9999 minutes.

Phase 4 - Exhaustion

Composed by releasing the pressure from the chamber and cooling the fluids that will be sent to the drainage system. No programming parameters.

Phase 5 - Drying

With a new vacuum phase to eliminate moisture from the chamber and dry the load. The parameters for this phase are:

- Drying pulses with a programming range between 0 and 9999.
- Vacuum level that assumes the value set for the initial vacuum pulses as standard.
- Drying time with programming range between 1 and 9999 minutes.

Phase 6 - Aeration

With the entry of filtered air through the hydrophobic filter, breaking the vacuum and equalizing the chamber pressure with the environment makes it possible to open the discharge door. No programmable parameters.

LIQUID CYCLES

Sterilization cycles with profiles for liquids in open bottles. With flexible programming the definition of parameters depends fundamentally on the characteristics of the liquid and the volume of the bottle where it is contained. These variables define the process temperature, the time and the characteristics of the initial phase of air removal.

The cycle allows the user to choose the load sensor, if one is available, for process control.

Phase 1 - Air removal and conditioning

For liquids in open flasks, air removal starts with a shallow vacuum pulse to prevent the liquid from

leaking out of the flasks and is complemented by a steam sweep, where the steam inlet and drain valve remain open, to let the steam drag the residual air in the chamber into the drain. Programmable parameters of this phase are:

- Single vacuum pulse with programmed depth in the range of 0,50 and 0,90 bar.
- Steam pulse with programming range between 1,0 and 2,0 bar.
- Scan, with programmable time between 01 and 9999 minutes.

Phase 2 - Load heating

With steam injection until reaching the sterilization temperature, in the chamber. Phase control through chamber temperature and pressure measurements.

Phase 3 - Exhibition

The time to reach temperature in the chamber basically depends on the volumes of liquid. When there is a core sensor, it is possible to determine that the start of the exposure time count will be defined by the measurements of that sensor. Without the sensor, the time to reach the temperature is defined in the process validation and added to the exposure time.

Programmable parameters of this phase are:

- Control Sensor Definition: between the sensor located in the drain and the load sensor in the chamber, when available.
- Sterilization temperature: with a programming range between 101 and 135°C.
- Exposure time: after reaching the temperature when using the control sensor, or exposure time as validated in the absence of the control sensor. With programming range between 01 and 9999 minutes.
- Cycle control by FO: with the definition of the desired lethality equivalent to an exposure to 121°C. With programming range between 01 and 9999.

Phase 4 - Cooling

With the slow release of pressure from the chamber to prevent the liquid from boiling, until the temperature set for the end of the cycle. It is a programmable parameter of this phase:

End of cycle temperature: with programming range between 80 and 99°C

STERILIZATION CYCLES AVAILABLE

Test cycles:

Cycle - Bowie & Dick

To check the air removal phase. With 4 vacuum pulses up to 0,20 bar, vapor pressure of 1,15 bar and a waiting time of 30 seconds, configurable in the available programming ranges.

Temperature of 134°C, exposure time of 210 seconds and drying time of 01 minutes, according to specific norm, without possibility of alteration.

Cycle - "Leak test"

To check the stability of the chamber. Non-programmable cycle, with factory fixed parameters. Vacuum pulse up to 0,20 bar, interval for accommodating the pressure of 05 minutes and monitoring the pressure from the stabilized level for 10 minutes. The test is considered "accepted", when the pressure increase is less than 0,013 bar.

Density profile cycles

Cycle - Textiles

Phases of air removal and conditioning, with the parameters defined for the B&D test. Temperature 134°C, exposure time of 04 minutes and drying time of 15 minutes, configurable in the available ranges of each parameter.

Cycle - Thermosensitive

Phases of air removal and conditioning, with the parameters defined for the B&D test. Temperature 121°C, exposure time of 20 minutes and drying time of 20 minutes, configurable in the available ranges of each parameter.

Cycle - Fast (flash)

For unpacked or soaked material. Phases of air removal and conditioning, with the parameters defined for the B&D test. Temperature 134°C, exposure time of 04 minutes and drying time of 05 minutes, configurable in the available ranges of each parameter.

Cycle - Instrumental

Phases of air removal and conditioning, with the parameters defined for the B&D test. Temperature 134°C, exposure time of 04 minutes and drying time of 30 minutes, configurable in the available ranges of each parameter.

Cycle - User A

It allows the configuration of the removal phase of independent values defined for Bowie & Dick. Temperature of 134°C, exposure time of 4 minutes and drying time of 30 minutes. All parameters, including those of the air removal phase, configurable in the permitted ranges for each parameter.

Cycles profile for liquids:

Cycle - Liquids

Single and initial vacuum pulse up to 0.50 bar; scanning pressure of 1,15 bar; scan time of 10 min; sterilization temperature of 121°C; time of 30 minutes and end of cycle temperature of 95°C. All configurable parameters in the ranges available for liquid cycles.

Cycle - User B

Configurable cycle for liquid profiling.

Lactate cycle

The parameters indicated for the lactate cycle are: single and initial pulse of 0,50 bar; 1,15 bar sweep pulse; 10 minute scan time; sterilization temperature 110°C, time of 10 minutes and end of cycle temperature of 95°C. All parameters configurable in the ranges available for liquid cycles.

OTHER INSTRUMENTS

Manovacuumeter for indicating the vacuum or pressure level of the inner chamber. With DN 63 mm bourdon tube, stainless steel system mechanism with a scale of 0 to 6 bar absolute, calibrated and filled with glycerin for stabilization.

Pressure gauges for pressure indication, one for an external chamber and a second for the steam generator. With DN 63 mm bourdon tube, stainless steel system mechanism with a scale of 0 to 6 bar absolute, calibrated and filled with glycerin for stabilization.

Thermal printer 40 column thermal printer, high resistance thermoplastic body and thermal printing technology, paper roll protected by cover, serial interface and installed on the front pa-

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The printed data guarantees a complete sterilization process with information such as:

- Name of the establishment; sterilizer identification; Dice; cycle start time; cycle number; charge code.
- Program number and name with parameters reached in the process. Cycle control information by FO, when applicable.
- Information about the phases performed, with time recording in minutes and seconds and data corresponding to the temperature of the control sensor, temperature of the verification sensor, pressure of the chamber transducer and available temperature of the test load sensor.
- No final information on the cycle and time of the process. When controlled by "FO", the corresponding "F" value will be printed.
- Failure to execute the cycle.

Matrix printer 40 columns composed of thermoplastic body of high strength and printing paper and the head coil protective cover with an impact printing mechanism quick eight needles. With serial interface and installed in the sterilizer front panel.

Alternatives to Printers

- Thermal printer
- Matrix printer

Process Control Valves with independent pneumatic drive, built in brass. Controlled by compressed air solenoid valves, in the steam supply lines of the two chambers, the vacuum and exhaust lines, and the air inlet at the end of the process.

Other controls such as pressurization and vacuum of port channels and water inlet for pumps through solenoid valves.

Control Valve Alternatives Pneumatic Valves

- Bronze alloys

OPTIONAL COMPONENTS STEAM GENERATION SYSTEM

Electric steam generator with automatic water supply, through a centrifugal pump with nickel chromium treatment. The generator body is constructed of AISI 316L stainless steel and armored heating element in stainless steel AISI 316.

Control of water level by stainless steel float system, for detection of minimum and maximum level of water in the generator.

Safety valve with calibration certificate, in AISI 316 stainless steel and threaded connections, sealed and with function check mechanism, protects the generator set. Built

according to current regulations.

Steam switch

System for steam power switching, electronically actuated through the sterilizer control interface. To replace the steam generator power to an embedded device to a central source of supply. Threaded connection stainless steel or brass as the material of choice for pipes of the equipment.

Alternatives steam switch

- With steam switch
- Without steam switch

Generator feed pump

The steam generator capacity influences the sterilization cycle time. To meet different demands and the availability of electrical power in the installation locations are different powers available to the same model of equipment.

Power steam generator

Electrical generators

- With 18 kW power.
For 100 l sterilizers.
- With a power of 24 kW.
For 210, 300 and 500 l sterilizers.
- With a power of 36 kW.
For 210, 300 and 500 l sterilizers.
- With a power of 44 kW.
For 210, 300 and 500 l sterilizers.

COMPRESSED AIR SUPPLY

The machine requires compressed air to drive the sealing ports and valves. For the supply of compressed air is available air compressor mounted outside of the machine structure. With the reservoir 50 or 100 liters and capable of supplying compressed air to 8,0 kgf / cm² pressure.

Alternatives for compressed air

- Compressed air network
- With air Compressor

REVERSE OSMOSIS

Optional system for treating feed water for steam generator, provided with three filter stages - with a degree of retention of 5 and 1 micron, and activated carbon - pressurizing pump, reverse osmosis membranes, control and treated water tank.

Alternatives to Reverse Osmosis

- Without osmosis system

- With osmosis system to 20 l/h for sterilizing 100-210 l

- With osmosis system of 60 l / h for sterilizing 300-500 l

SAFETY

The equipment has the following safety items:

- In the absence of electricity, the steam inlet is closed;
- When there is excess pressure, the vapor passage is closed;
- Safety valve standardized and calibrated at 3.0 kgf/cm² with cleaning device and function check;
- Start of the cycle only in the presence of steam pressure in the generator or in the supply line;
- Automatic door lift system with calibrated working pressure, to protect the operator;
- Door with anti-crushing system that interrupts the movement of the door if it is obstructed;
- Impossibility of opening the doors after the beginning of the cycle and in the presence of vapor pressure in the internal chamber;
- The cycle starts only with closed doors and pressurized linings;
Impossibility of opening the two doors simultaneously by the operator (in the case of barrier sterilizers);
- Abortion of the cycle when the chamber temperature in the sterilization phase is below that defined for the process;
- Safety thermostat to protect electrical resistances;
- Emergency routine that, in the absence of electricity, keeps the doors from locking and the pressure release from the internal chamber;
- Alarm system for water, steam and compressed air supply failure;
- Emergency button on the equipment panel to switch off the sterilizer and depressurize the sterilization chamber;
- Mechanical lock that prevents the door from opening when the sterilizer is in operation.

Note

The total width of the unpacked equipment is indicated in the engineering data table. To facilitate the movement of the equipment inside the building, the front can be disassembled, reducing the width. It is important to check the route that will be taken by the sterilizer inside the building to check for bottlenecks. The release of these points is the buyer's responsibility;

The gauges of the connections and supply lines indicated in technical data refer to the connection points with the sterilizer. The pipes that take these supplies to the installation sites must be calculated by third parties and the dimensioning depends on the distance, supply pressures and necessary flows.

ENGINEERING DATA

DRAIN

The drain line must be independent of the sewer building and connected to an external junction box. The connection of the equipment to the building sewer line is only indicated when the equipment has heat exchanger to cool the effluents released in the process. Piping made of heat resistant material.

WATER SUPPLY

The sterilizers have water inlet to the steam generator and the vacuum system. In generator for water inlet can be connected to a reverse osmosis system to treat water used for steam generation.

Water supply for Vacuum Pump

- Maximum hardness 250 ppm CaCO₃.
- Ideal temperature lower than 15°C. Maximum allowable 25°C.
- pH 5 to 8.
- Pressure between 1.2 and 1.5 bar absolute.

Water supply for reverse osmosis

- Hardness less than 60 ppm CaCO₃.
- Conductivity less than 500 µS/cm.
- Water temperature less than 25°C
- pH 5 to 8.
- Chloride <60 mg/l.
- Dynamic pressure: min 1 bar and a maximum of 4 bar.

Water supply for steam generator

- Hardness less than 10 ppm CaCO₃.
- Conductivity less than 50 µS/cm
- Water temperature between 5 and 35°C
- pH 5 to 8.
- Chloride <3 mg/l.
- Dynamic Pressure: Minimum 1 bar absolute.

NOTE: The required flow rates vary with the configuration and model of equipment. See the Installation Guide for more information.

STEAM SUPPLY

The direct steam supply to the autoclave chamber is not suitable solution for meeting the requirements of Good Practices for material sterilization processes.

Direct Steam supply

- Quality suggestion as EN 285 and ISO 17665-2.
- Pressure 3,5 to 3,7 bar absolute.
- Steam title above 0,95.

NOTE: The Baumer is not responsible for stains on instruments and tissue and corrosion on instruments and pieces of equipment in case of direct steam supply outside the standards suggested by the standard above.

The best alternative is the provision of process steam for a steam Generator with steam coil.

Steam supply for Steam Generator with steam coil

- Pressure 7-9 bar absolute.
- Steam title above 0,95.

COMPRESSED AIR SUPPLY

Compressed air free of impurities and moisture and supply pressure of 6 to 7 bar. Push-in connector 6 mm.


HI SPEED II line has an internal electrical panel with three-phase and earth connection with a maximum impedance of 10 Ohms and protections for each load element per circuit breaker.

NOTE: Refer to the Installation Guide to identify

ELECTRICAL SUPPLY

EQUIPMENT CAPACITY AND DIMENSIONING

| Capacity (L) | 100 | 210 | 300 | 500 |
|--|--------------|-----------------|---------------|-----------------|
| model 1P | B0111-100 | B0111-210 | B0111-300 | B0111-500 |
| model 2P | B0111-100P | B0111-210P | B0111-300P | B0111-500P |
| Dim. Chamber (cm) | 40,8x40,8x63 | 40,8x40,8x128,5 | 60,5x60,5x88 | 60,5x60,5x128,5 |
| Cap. Baskets ISO | 2 | 4 | 6 | 9 |
| A - Width (mm) | 920 | 920 | 1195 | 1195 |
| B - Depth (mm) | 985 | 1585 | 1180 | 1585 |
| C - Height (mm) | 1710 | 1710 | 1880 | 1880 |
| BF - Free space for loading and unloading (mm) | 1750 | 2000 | 1750 | 2000 |
| Packaging (cm) | 108x119x200 | 108x178,5x200 | 135,5x138x217 | 135,5x178x217 |
| Net Weight (kg) | 350 | 830 | 850 | 875 |
| Weight (kg) | 500 | 960 | 990 | 1075 |

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