

Data Sheet

Reusable Collection Jar



Product Name	Reusable Collection Jar
Manufacturer	Gentec (Shanghai) Corporation
Part No.	107100150
Model No.	SJ-1500
Material	Silicone, PC, Steel 304
Straight Connector	Hose Barb(ϕ 6-8)
Intended use	The suction jar is used to collect pus, blood clots, phlegm, and excess fluids by connecting to a patient on one end, and vacuum source (such as suction regulator) on the other.
Storage conditions	Storage Temperature: -20°C~ 60°C Storage Relative Humidity: < 90%
Weight (kg)	0.6
Dimensions	60x45.5x26cm
Qty.(pcs/box)	12

Reusable Collection Jar



Product Name	Reusable Collection Jar
Manufacturer	Gentec (Shanghai) Corporation
Part No.	107100253
Model No.	SJ-2500
Material	Silicone, PC, Steel 304
Straight Connector	Hose Barb(φ6-8)
Intended use	The suction jar is used to collect pus, blood clots, phlegm, and excess fluids by connecting to a patient on one end, and vacuum source (such as suction regulator) on the other.
Storage conditions	Storage Temperature: -20°C~ 60°C Storage Relative Humidity: < 90%
Weight (kg)	0.7
Dimensions	67×53.4×33.5cm
Qty.(pcs/box)	12