

NEBULIZER KIT



637-003 Series

Solutions for Life




Features

- NebuSmart is All-in-one Nebulizer :
A: Breath-enhance B. Breath-closed C Breath-actuated
For treatment of upper, middle, and lower airways
- Anti-spill design. Can be tilted 90° vertically without spilling
- Clear calibration marks for easy delivery of medicines
- Single patient use to avoid cross infection
- Easy and fast to deliver the drugs during nebulization (figure 1)
- Breath- closed is used for beathing circuit with T-piece (figure 2)




A. Breath-enhanced

Faster and higher output for the treatment of upper airway



B. Breath-closed

closed system with smallest partical size for the treatment of lower airway



C. Breath-actuated

one way valve to increase the efficiency of treatment

Ordering Information



Kid Set



Mouth Piece Set



Adult Set



Pacifier Set

Part Number	Model Number	Description							Specifications
		Oxygen Tube	Jar	Seal	One Way Valve	Angle Mouth Piece	Mask	T-piece	
A63202320	637-003-S08	✓	✓	✓	✓	-	S	-	Kid Set
A63202330	637-003-M08	✓	✓	✓	✓	✓	-	-	Mouth Piece Set
A63202340	637-003-L08	✓	✓	✓	✓	-	L	-	Adult Set
A63202350	637-003-P08	✓	✓	✓	✓	-	-	-	Pacifier Set
A63202360	637-003-A08	✓	✓	✓	✓	-	✓	✓	Universal Set