

Soft, transparent and with adjustable nose clip, GENTEC® oxygen masks come with kink resistant oxygen tube and reservoir bag (optional).



Name	Model	Specifications
Oxygen Mask (Medium Concentration)	642-11C	Child: 5mmID tube, length: 2m
	642-11S	Adult S: 5mmID tube, length: 2m
	642-11M	Adult M: 5mmID tube, length: 2m
	642-11L	Adult L: 5mmID tube, length: 2m



Name	Model	Specifications
Oxygen Mask (High Concentration)	642-21C	Child: with reservoir bag, single patient use
	642-21S	Adult S: with reservoir bag, single patient use
	642-21M	Adult M: with reservoir bag, single patient use
	642-21L	Adult L: with reservoir bag, single patient use



Name	Model	Specifications
Oxygen Mask (Two color variable Concentration)	642-31C	Green diluter for 24%, 26%, 28%, 30% With diluter for 35%, 40%, 50%
	642-31S	
	642-31M	
	642-31L	



Name	Model	Specifications
Oxygen Mask (Seven color variable Concentration)	642-32C	Blue 24%, yellow 28%, white 31%, green 35%, pink 40%, orange 50%, purple 60%
	642-32S	
	642-32M	
	642-32L	