

Electrosurgical Unit

ERMCM 400 WATTS
ESU400

Power Supply:
110V, 50/60Hz or 220V, 50/60Hz

Dimension:
Product: 421(L) x 357(W) x 162(H) mm
Package: 520(L) x 445(W) x 305(H) mm

Weight:
Net Weight: 8.7KG
Gross Weight: 14.5KG




Standard Accessory



Optional




Note for above products:
Special length, tip and/or colour
can be custom ordered as per the demand



ERMCM-ESU400

High-frequency Electrocautery Unit

Easy-to-use buttons for simple operation
High power output 400W for more operation selection
Up to 7 work modes: Pure/Blend1/Blend2, Standard/Spray, Bipo1/Bipo2
Increment of 5W micro power control system to ensure optimal effect Minimal
sticking and carbonization, reduce the surgery damage
Fast trouble-shooting by default error alarm system
Auto-terminate of power output after 5 seconds to avoid misoperation Support
endoscopes connection and underwater cutting (TURP etc) OEM projects



Multiple Functions
Clear Color Indicators for Work
Mode & Alarm



Easy-to-Use
Convenient Operation by
Finger-suit Buttons



Error Code
Simple Error Code for Fast
Trouble-shooting



High Safety
Rem(Return Electrode
Monitoring) system to avoid
any burnt or damage to the
surgeon and patient.



Wide Application
It can be connected to
endoscopes for more clinical
applications.

Note for above products:
Special length, tip and/or colour
can be custom ordered as per the demand

