


# REUSABLE BREATHING CIRCUIT



## Features

- One piece extrusion made by high tension medical grade silicone
- Corrugated smoothbore tubing
- Smooth inner bore allows for less resistance and reduced turbulence
- Ultra-flexible with low circuit compliance and resistance free
- Kink resistant spiral outer casing
- Crush resistant design allows for flow even when twisted through 180°
- Standard tube sizes ID: 19F-Adult, 15F-Pediatric, 10F-neonatal
- High flexibility with low circuit compliance and resistance free
- Additional movement range for patient through use of a longer Catheter mount
- Non adhesive cuff connection increases durability and reliability
- ISO standard connectors
- Water trap can be completely disassembled for cleaning
- All parts and components are autoclavable (up to 134°C)
- Wye piece can be connected with GENTEC NebuCare in-line electrical nebulizer



663 Series  
Anesthesia, Breathing Circuit for Adult

## Dimensions

Type	A Water trap (PSF)	B (Smoothbore Tube)			C Swivel Elbow with Suction Port (PSF)	D Catheter Mount	E Piece with Temperature & Pressure Port
		B1	B2	B3			
Adult	40 ml 22 M / 22 M	22 F, 22 F 600mm(24"), 19 ID	22 F, 22 F 450mm(18"), 19 ID	22 F, 22 F 1200mm(48"), 19 ID	22 M / 15 F 15 M	22 F / 15 F 180mm(7")	22 M / 15 F 22 M
Pediatric	40 ml 22 M / 22 M	22 F, 22 F 600mm(24"), 15 ID	22 F, 22 F 450mm(18"), 15 ID	22 F, 22 F 1200mm(48"), 15 ID	22 M / 15 F 15 M	22 F / 15 F 180mm(7")	22 M / 15 F 22 M
Neonate	40 ml 15 M / 15 M	10 F, 15 F 600mm(24"), 15 ID	10 F, 15 F 450mm(18"), 15 ID	10 F, 10 F 1300mm(48"), 15 ID	15 M / 15 M	10 F / 15 F 350mm(14")	10 M / 15 M