

Data Sheet

Face Mask



Product Name	NIV Full Face Mask (Forehead Rest),L
Manufacturer	Gentec (Shanghai) Corporation
Part No.	A62200110
Model No.	641-008-NL
Material	Silicone, PC
Straight Connector	22OD
Intended use	This product is connected to the positive pressure ventilation ventilator through the breathing line, and the cover body is fixed on the user's nose, mouth or face, so that the positive pressure air flows through the nose and/or mouth into the airway, which serves as the interface for the patient to receive treatment gas.
Storage conditions	Storage Temperature: -20°C~ 50°C Storage Relative Humidity: < 80%
Weight (kg)	0.2
Dimensions	56x52x27cm
Qty.(pcs/box)	50

Face Mask



Product Name	NIV Full Face Mask (Forehead Rest),M
Manufacturer	Gentec (Shanghai) Corporation
Part No.	A62200070
Model No.	41-008-NM
Material	Silicone, PC
Straight Connector	22OD
Intended use	This product is connected to the positive pressure ventilation ventilator through the breathing line, and the cover body is fixed on the user's nose, mouth or face, so that the positive pressure air flows through the nose and/or mouth into the airway, which serves as the interface for the patient to receive treatment gas.
Storage conditions	Storage Temperature: -20°C~ 50°C Storage Relative Humidity: < 80%
Weight (kg)	0.2
Dimensions	56x52x27cm
Qty.(pcs/box)	50